QUICK START GUIDE

SMART SCALE
8 in 1

Model: HSCALE2

2	English GB
4	Polish PL
7	German DE
10	Greek EL
14	Czech CZ
19	Hungarian HU
22	Italian IT
25	Lithuanian LT
31	Romanian RO
34	Slovenian SI
37	Slovak SK
39	Ukrainian UA

Device Overview - FIG 1


FIG 2 - PL




FIG 2 - GB

Device Overview - FIG 1

1. LCD display
2. Electrode
3. Battery Compartment

How to use scale with a mobile APP? – EX. II

1. Download „OK OK – International“ APP
2. If you do not have a personalised account yet, please create one (1A). If you already have an account, please log in (1B).
3. In order to get optimal results, please fill in personalised information about yourself (2A/3).
4. Turn on Bluetooth in device settings. If you use Android newer than 5.1, you can meet the request to allow location settings. Pair the scale with your device (2B/4). The scale will get paired, if the user stands on it barefoot and there will be there will be weight result shown on LCD screen. Make sure the LCD screen and APP show the same result and confirm.
5. The scale needs a few seconds to analyse all parameters, do not step down the scale before they are loaded (2D/5). Parameters will only be loaded if you stand on the scale barefoot.
6. Main screen of the APP will show its user latest result, you can explore full body analysis by clicking on the weight result number (2C/5). The user can also choose older results from the list (6).

BIA (Bio-Electrical Impedance Analysis)

1. GOCLEVER SCALE lets a small amount of weak (less than 1mA) current flow through the human body, so as to detect the bioimpedance and estimate body fat. However, please be aware that anyone with an wearable or implantable medical electronic instrument, such as a pacemaker, must avoid using this device.

Disclaimer

Please be aware that this device is designed for adults' self-measuring and self-monitoring body fat level. Any information provided by this device is in no way meant to treat, cure or prevent any disease or illness from happening. This device should not be used by anyone who is acutely or chronically ill, suffering from a disease or taking medications that affect your water levels.

DECLARATION OF CONFORMITY

GOCLEVER Ltd. hereby declares under our sole responsibility that the product GOCLEVER SMART SCALE 8in1 (HSCALE2) is in conformity with the essential requirements of the following Directives is declared:

EMC DIRECTIVE 2004 / 108 / EC


R&TTE DIRECTIVE 1999 / 5 / EC

RoHS DIRECTIVE 2011/65/EU

The full document with detailed information is available on our website: www.goclever.com

Correct Disposal of This product (Waste Electrical & Electronic Equipment). (Applicable in the European Union and other European countries with separate collection systems) This marking indicates that at the end of its working life this product should not be disposed of with other household waste in the EU. To prevent possible harm to the environment or human health from uncontrolled waste disposal, please separate this from other types of wastes and recycle it responsibly to promote the sustainable reuse of material resources. To dispose of your used device, please use the return and collection systems available in your area or contact the retailer where the product was purchased. They can take this product for safe environmental recycling. Correct Disposal of This product (Waste Electrical & Electronic Equipment).

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